

(R,S)-Cotinine-D3 N-Oxide

Chemical Properties


CAS No. : 1215561-37-7

Formula: C₁₀H₁₂N₂O₂

Molecular Weight: 195.23

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(R,S)-Cotinine-D3 N-Oxide is the deuterated form of (R,S)-Cotinine N-Oxide.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.1222 mL	25.6108 mL	51.2216 mL
5 mM	1.0244 mL	5.1222 mL	10.2443 mL
10 mM	0.5122 mL	2.5611 mL	5.1222 mL
50 mM	0.1024 mL	0.5122 mL	1.0244 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481