

## Pyloricidin C

## Chemical Properties


CAS No. : 219823-30-0

Formula: C<sub>21</sub>H<sub>33</sub>N<sub>3</sub>O<sub>8</sub>

Molecular Weight: 455.50

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Pyloricidin C is an antibiotic that targets Helicobacter pylori.
Targets(IC50)	Antibacterial,Antibiotic

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1954 mL	10.9769 mL	21.9539 mL
5 mM	0.4391 mL	2.1954 mL	4.3908 mL
10 mM	0.2195 mL	1.0977 mL	2.1954 mL
50 mM	0.0439 mL	0.2195 mL	0.4391 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481