

PRMT5-IN-52

## Chemical Properties


CAS No. : 2783961-64-6

Formula: C19H16N4O

Molecular Weight: 316.36

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	PRMT5-IN-52 (Compound 17) is a non-nucleoside inhibitor of protein arginine methyltransferase 5 (PRMT5), demonstrating an inhibition rate of 20.2% at a concentration of 10 $\mu$ M. PRMT5-IN-52 exhibits potent antitumor activity and holds potential for research in cancers such as lung, prostate, and colorectal cancer.
Targets(IC50)	Histone Methyltransferase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.161 mL	15.8048 mL	31.6096 mL
5 mM	0.6322 mL	3.161 mL	6.3219 mL
10 mM	0.3161 mL	1.5805 mL	3.161 mL
50 mM	0.0632 mL	0.3161 mL	0.6322 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481