


Kv7.2 activator-1

Chemical Properties

CAS No. :	2981435-59-8
Formula:	C19H23FN4O2
Molecular Weight:	358.41
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Kv7.2 activator-1 (Compound 8) is an activator of Kv7.2/3 channels with an EC50 of 178 nM. It has a relatively long half-life of 120 minutes in human liver microsomes (HLM).
Targets(IC50)	Potassium Channel

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7901 mL	13.9505 mL	27.901 mL
5 mM	0.558 mL	2.7901 mL	5.5802 mL
10 mM	0.279 mL	1.3951 mL	2.7901 mL
50 mM	0.0558 mL	0.279 mL	0.558 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481