


Piptamine

Chemical Properties

CAS No. :	5001-21-8
Formula:	C ₂₀ H ₃₅ N ₃
Molecular Weight:	317.51
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Piptamine exhibits a broad antibacterial spectrum, with inhibitory MIC ($\mu\text{g}/\text{mL}$) values against various microorganisms as follows: Staphylococcus aureus (0.78-6.25), Enterococcus faecalis (1.56), Bacillus subtilis (1.00), Escherichia coli (12.5), and Candida albicans (6.25).
Targets(IC50)	Antibacterial

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1495 mL	15.7475 mL	31.4951 mL
5 mM	0.6299 mL	3.1495 mL	6.299 mL
10 mM	0.315 mL	1.5748 mL	3.1495 mL
50 mM	0.063 mL	0.315 mL	0.6299 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481