


RNase L ligand 3

Chemical Properties

CAS No. :	3036649-51-8
Formula:	C27H27N3OS
Molecular Weight:	441.59
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	RNase L ligand 3 is an RNase L ligand employed in the synthesis of F3-PEG8-RiboTAC.
Targets(IC50)	DNA/RNA Synthesis,PROTACs

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2645 mL	11.3227 mL	22.6454 mL
5 mM	0.4529 mL	2.2645 mL	4.5291 mL
10 mM	0.2265 mL	1.1323 mL	2.2645 mL
50 mM	0.0453 mL	0.2265 mL	0.4529 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481