


3-Methyldodecanoyl-CoA

Chemical Properties

CAS No. :	711027-67-7
Formula:	C34H60N7O17P3S
Molecular Weight:	963.86
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-Methyldodecanoyl-CoA (3-Methyldodecanoyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0375 mL	5.1875 mL	10.375 mL
5 mM	0.2075 mL	1.0375 mL	2.075 mL
10 mM	0.1037 mL	0.5187 mL	1.0375 mL
50 mM	0.0207 mL	0.1037 mL	0.2075 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481