

3-Oxononanoyl-CoA

Chemical Properties


CAS No. : 1206185-44-5

Formula: C₃₀H₅₀N₇O₁₈P₃S

Molecular Weight: 921.74

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-Oxononanoyl-CoA (3-Oxononanoyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0849 mL	5.4245 mL	10.849 mL
5 mM	0.217 mL	1.0849 mL	2.1698 mL
10 mM	0.1085 mL	0.5425 mL	1.0849 mL
50 mM	0.0217 mL	0.1085 mL	0.217 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481