


IRAK4 ligand-14

Chemical Properties

CAS No. :	3065496-26-3
Formula:	C ₂₆ H ₃₀ N ₄ O ₄
Molecular Weight:	462.54
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	IRAK4ligand-14 is a ligand for the IRAK4 protein. It is useful for synthesizing PROTACs, such as [APH02174].
Targets(IC50)	Ligands for Target Protein for PROTAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.162 mL	10.8099 mL	21.6198 mL
5 mM	0.4324 mL	2.162 mL	4.324 mL
10 mM	0.2162 mL	1.081 mL	2.162 mL
50 mM	0.0432 mL	0.2162 mL	0.4324 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481