


## 8-Oxo-dGTP

## Chemical Properties

CAS No. :	139307-94-1
Formula:	C <sub>10</sub> H <sub>16</sub> N <sub>5</sub> O <sub>14</sub> P <sub>3</sub>
Molecular Weight:	523.18
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	8-Oxo-dGTP is an oxidized guanine nucleotide and serves as a critical biomarker of oxidative stress.
Targets(IC50)	ROS

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9114 mL	9.5569 mL	19.1139 mL
5 mM	0.3823 mL	1.9114 mL	3.8228 mL
10 mM	0.1911 mL	0.9557 mL	1.9114 mL
50 mM	0.0382 mL	0.1911 mL	0.3823 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481