

Pomalidomide-NH-C13-OH

Chemical Properties

CAS No. : 2716120-15-7


Formula: C₂₆H₃₇N₃O₅

Molecular Weight: 471.59

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pomalidomide-NH-C13-OH is a synthetic E3 ligase ligand-linker conjugate (E3 ligase Ligand-linker conjugate) that is utilized in the synthesis of PROTAC PI3K α degrader-1.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1205 mL	10.6024 mL	21.2049 mL
5 mM	0.4241 mL	2.1205 mL	4.241 mL
10 mM	0.212 mL	1.0602 mL	2.1205 mL
50 mM	0.0424 mL	0.212 mL	0.4241 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481