


## Padoprazanum fumarate

## Chemical Properties

CAS No. :	2938973-98-7
Formula:	C <sub>23</sub> H <sub>24</sub> FN <sub>3</sub> O <sub>8</sub> S
Molecular Weight:	521.52
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Padoprazanum (Padoprazan) fumarate is a proton pump (Proton Pump) inhibitor.
Targets(IC50)	Proton pump

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9175 mL	9.5874 mL	19.1747 mL
5 mM	0.3835 mL	1.9175 mL	3.8349 mL
10 mM	0.1917 mL	0.9587 mL	1.9175 mL
50 mM	0.0383 mL	0.1917 mL	0.3835 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481