


GLP-1R agonist 31

Chemical Properties

CAS No. :	2922542-61-6
Formula:	C ₃₂ H ₃₁ Cl ₂ N ₅ O ₅
Molecular Weight:	636.53
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	GLP-1R agonist 31 (Compound 1) is an amorphous glucagon-like peptide-1 receptor (GLP-1R) agonist.
Targets(IC50)	Glucagon Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.571 mL	7.8551 mL	15.7102 mL
5 mM	0.3142 mL	1.571 mL	3.142 mL
10 mM	0.1571 mL	0.7855 mL	1.571 mL
50 mM	0.0314 mL	0.1571 mL	0.3142 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481