

ALD1910

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ALD1910 is a humanized monoclonal antibody targeting PACAP38 and PACAP27. It inhibits PACAP signaling through pituitary adenylate cyclase-activating peptide type I receptor (PAC1-R), vasoactive intestinal peptide receptor 1 (VPAC1-R), and VPAC2-R. ALD1910 recognizes a non-linear epitope within PACAP and prevents its binding to cell surfaces. It is useful for studying PACAP-mediated migraines.
Targets(IC50)	PACAP

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481