

ADI-65533

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	ADI-65533 is a human IgG1 monoclonal antibody that targets Envelopment polyprotein, M polyprotein. It is useful for research on hantavirus infections. The recommended isotype control is Human IgG1 kappa, Isotype Control.
Targets(IC50)	Influenza Virus

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481