

FAZ-053

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	FAZ-053 is a human monoclonal antibody (mAb) that targets B7-H1/PD-L1/CD274. It inhibits the interaction of PD-L1 with PD-1 and B7-1 on monocytes, dendritic cells, and B cells. Additionally, FAZ-053 enhances the production of interleukin-2 and is applicable in research studies on advanced alveolar soft part sarcoma (ASPS), chordoma, and triple-negative breast cancer.
Targets(IC50)	PD-1/PD-L1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481