

IBI37G5

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	IBI37G5 is a human monoclonal antibody that targets TNFRSF18/GITR/CD357. It demonstrates antitumor activity in MC38 and B16F10 mouse tumor models. IBI37G5 is applicable in cancer immunology research.
Targets(IC50)	TNF

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481