

ZEB85

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at low temperature

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ZEB85 is a human monoclonal antibody (mAb) that targets TrkB. It activates TrkB and its downstream signaling cascade, enhancing neuronal activity. ZEB85 can prevent β -amyloid toxicity in cultured hippocampal neurons and is useful for Alzheimer's disease (AD) research.
Targets(IC50)	Trk receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481