

BIIB022

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	BIIB022 is a human IgG4 monoclonal antibody (mAb) targeting IGF1R/CD221. It inhibits the phosphorylation of IGF-1 and IGF-2-induced IGF-1R and its downstream substrates. BIIB022 is applicable for research in lung cancer, pancreatic cancer, and colon cancer. The recommended isotype control is Human IgG4 kappa, Isotype Control.
Targets(IC50)	IGF-1R

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481