

IDB-002

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	IDB-002 is a human monoclonal antibody (mAb) targeting the F/FusionglycoproteinF0. It is suitable for research related to respiratory syncytial virus infections. The recommended isotype control for this study is Human IgG1 kappa, Isotype Control.
Targets(IC50)	RSV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481