

ASKG-712

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ASKG-712 (AM712) is a recombinant humanized monoclonal antibody (mAb) combined with an Ang-2 antagonist peptide fusion protein. It is utilized in research focused on wet age-related macular degeneration and diabetic macular edema.
Targets(IC50)	Tie-2,VEGFR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481