

VRDN-002

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	VRDN-002 is a humanized monoclonal antibody (mAb) targeting IGF1R/CD221. It is designed to reduce inflammation and exophthalmos associated with thyroid eye disease (TED). Additionally, VPI-2690B is utilized in research related to Graves' ophthalmopathy.
Targets(IC50)	IGF-1R

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481