

ADI-42898

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	ADI-42898 is a human IgG1 monoclonal antibody targeting Envelopment polyprotein and M polyprotein. It binds to the Gn/Gc complex, preventing the entry of seven Hantavirus strains associated with human immunodeficiency syndrome (HCPS) and hemorrhagic fever with renal syndrome (HFRS) into cells. Recommended isotype control: HumanIgG1kappa, Isotype Control.
Targets(IC50)	Influenza Virus

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481