


## 3-Ketostearoyl-CoA

## Chemical Properties

CAS No. :	86370-20-9
Formula:	C39H68N7O18P3S
Molecular Weight:	1047.98
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-Ketostearoyl-CoA (3-Ketostearoyl-Coenzyme A) is a type of 3-oxo-fatty acyl-CoA.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9542 mL	4.7711 mL	9.5422 mL
5 mM	0.1908 mL	0.9542 mL	1.9084 mL
10 mM	0.0954 mL	0.4771 mL	0.9542 mL
50 mM	0.0191 mL	0.0954 mL	0.1908 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481