

SEA-TGT

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	SEA-TGT is a human monoclonal antibody (mAb) targeting TIGIT. It enhances binding to activating FcγRIIIA while reducing binding to the inhibitory FcγRIIb. SEA-TGT is applicable for research in advanced solid tumors and lymphomas.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481