


Acetyl-L-leucyl-L-argininal

Chemical Properties

CAS No. : 35039-75-9
 Formula: C₁₄H₂₇N₅O₃
 Molecular Weight: 313.40
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Acetyl-L-leucyl-L-argininal inhibits the activity of peptidylaminopeptidase III.
Targets(IC50)	Proteasome

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1908 mL	15.9541 mL	31.9081 mL
5 mM	0.6382 mL	3.1908 mL	6.3816 mL
10 mM	0.3191 mL	1.5954 mL	3.1908 mL
50 mM	0.0638 mL	0.3191 mL	0.6382 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481