

VPI-2690B

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	VPI-2690B is a human-derived monoclonal antibody (mAb) that targets Integrin $\alpha V\beta 3$. It is capable of modulating TGF- β signaling and can be utilized in research related to diabetic nephropathy.
Targets(IC50)	Integrin,TGF-beta/Smad

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481