

BIIB036

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | BIIB036 is a human-derived IgG1 monoclonal antibody (mAb) that targets TNFRSF12A/Fn14. It effectively inhibits the growth of patient-derived primary colorectal tumors in SCID mice. |
| Targets(IC50) | TNF |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481