

PE0116

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	PE0116 is a fully humanized CD137 agonistic monoclonal antibody, generated from immunized Harbor H2L2 transgenic mice. Acting as a ligand blocker, PE0116 activates NF-κB signaling, which significantly enhances T cell proliferation and increases cytokine secretion in vitro. In the MC38 tumor model, PE0116 demonstrates strong antitumor activity.
Targets(IC50)	NF-κB

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481