

BI-836858

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	BI-836858 is a fully humanized IgG1 anti-CD33 monoclonal antibody. It is modified at the Fc fragment to enhance binding with FcγRIIIa (CD16). BI-836858 exhibits antibody-dependent cellular cytotoxicity (ADCC) activity and is utilized in research related to acute myeloid leukemia (AML).
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481