

MEDI-1912

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	MEDI-1912 is a human monoclonal antibody targeting NGF/bNGF. It inhibits signaling via the TrkA and p75 receptors and can be utilized in chronic pain research.
Targets(IC50)	Neurokinin receptor,Trk receptor

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481