

NH₂-PEG2-COOMe

Chemical Properties


CAS No. : 2445792-36-7

Formula: C₈H₁₇NO₄

Molecular Weight: 191.23

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	NH ₂ -PEG2-COOMe is a biochemical reagent used in drug delivery research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2293 mL	26.1465 mL	52.2931 mL
5 mM	1.0459 mL	5.2293 mL	10.4586 mL
10 mM	0.5229 mL	2.6147 mL	5.2293 mL
50 mM	0.1046 mL	0.5229 mL	1.0459 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481