

## Anti-CD27 Antibody (M2177)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-CD27 Antibody (M2177) is a human-derived monoclonal antibody (mAb) that targets TNFRSF7/CD27. It can inhibit the binding of sCD70 to human CD27ECD-Fc and is applicable in anti-tumor immunology studies. Recommended isotype control: Human IgG1 kappa, Isotype Control.
Targets(IC50)	TNF

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481