

TRS-005 Antibody

Chemical Properties

CAS No. : 2851061-85-1

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	TRS-005 Antibody targets CD20-positive tumor cells and is suitable for synthesizing antibody-drug conjugates (ADC). It is also utilized in research on relapsed/refractory diffuse large B-cell lymphoma (R/R DLBCL).
Targets(IC50)	ADC Antibody,Immunology/Inflammation related

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481