

IMM20324

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	IMM20324 is a human-derived monoclonal antibody (mAb) targeting IL38. It inhibits the interaction of IL-38 with human IL1RAPL1-Fc and human IL-36R-Fc, with IC50 values of 1.267nM and 1.4667nM, respectively. In B16.F10 syngeneic mouse models, IMM20324 demonstrates antitumor activity and can be utilized in antitumor immunity research. The recommended isotype control is HumanIgG1kappa.
Targets(IC50)	Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481