

Sym004

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Sym004 is a 1:1 mixture of two recombinant human-mouse chimeric monoclonal antibodies (mAb), targeting non-overlapping EGFR epitopes (mAb992 and mAb1024). This compound facilitates the rapid internalization of EGFR and mediates subsequent degradation of internalized receptors through EGFR cross-linking. Sym004 can be utilized in research concerning metastatic colorectal cancer, head and neck cancer (H&N), and lung cancer.
Targets(IC50)	EGFR

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481