


5-Benzyloxytryptamine

Chemical Properties

CAS No. :	20776-45-8
Formula:	C17H18N2O
Molecular Weight:	266.34
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-Benzyloxytryptamine (5-BT) is a selective partial agonist of the 5-hydroxytryptamine 1D/1B (5-HT1D/1B) receptors, with an IC50 value of 40 nM for the bovine caudate nucleus 5-HT1D. It exhibits much lower affinity for other receptors (5-HT2 IC50 > 470 nM). 5-Benzyloxytryptamine inhibits adenylate cyclase, leading to reduced neurotransmitter release and downregulation of cAMP signaling. It holds potential for research into diseases associated with neurotransmitter imbalances.
Targets(IC50)	5-HT Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7546 mL	18.773 mL	37.546 mL
5 mM	0.7509 mL	3.7546 mL	7.5092 mL
10 mM	0.3755 mL	1.8773 mL	3.7546 mL
50 mM	0.0751 mL	0.3755 mL	0.7509 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481