

MIL62

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	MIL62 is a glycoengineered type II anti-CD20 monoclonal antibody with an Fc region N-glycan that is almost entirely afucosylated. It exhibits increased affinity for the FcγRIIIa receptor and possesses direct B-cell killing effects. MIL62 displays antibody-dependent cell-mediated cytotoxicity (ADCC) and is applicable for research in non-Hodgkin's lymphoma (NHL).
Targets(IC50)	Immunology/Inflammation related

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481