

NP137

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	NP137 is a human IgG1 monoclonal antibody (mAb) that targets Netrin-1. It inhibits the binding of netrin-1 to its receptor UNC5H2 with an IC50 of 0.5 nM. NP137 induces apoptosis by activating caspase 3 and is applicable for research in both solid tumors and hematological malignancies. Recommended isotype control: Human IgG1 kappa, Isotype Control.
Targets(IC50)	Apoptosis

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481