

BI-655064

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | BI-655064 is a humanized anti-CD40 antibody with two mutations in its crystallizable fragment (Fc) region, which prevent Fc-mediated antibody-dependent or complement-mediated cytotoxicity and platelet activation. It is utilized in the study of autoimmune diseases, including lupus nephritis and rheumatoid arthritis (RA). |
| Targets(IC50) | TNF |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481