

DS-1501

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at low temperature

Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DS-1501 is a human monoclonal antibody (mAb) that targets Siglec-15/CD33L3. It inhibits the formation of osteoclasts and prevents the decline of lumbar vertebrae bone mineral density (BMD) in ovariectomized (OVX) rats. DS-1501 is applicable in osteoporosis research. Recommended isotype control: HumanIgG1kappa [Isotype Control].
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481