

BI-836826

Chemical Properties

CAS No. : 2235420-23-0

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	BI-836826 is a chimeric IgG1 anti-CD37 monoclonal antibody with an engineered Fc fragment. It exhibits both antibody-dependent cellular cytotoxicity (ADCC) and direct pro-apoptotic activity. BI-836826 is applicable in research related to chronic lymphocytic leukemia.
Targets(IC50)	Apoptosis

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481