

RG-7258

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	RG-7258 is a humanized monoclonal antibody (mAb) that targets CRLF2. It inhibits dendritic cell activation and mast cell cytokine release, and reduces Th2 inflammation in non-human primate models of allergic pulmonary inflammation. RG-7258 is applicable in asthma research.
Targets(IC50)	Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481