


## 3-Oxodecanoyl-CoA

## Chemical Properties

CAS No. : 50411-91-1  
 Formula: C<sub>31</sub>H<sub>52</sub>N<sub>7</sub>O<sub>18</sub>P<sub>3</sub>S  
 Molecular Weight: 935.77  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-Oxodecanoyl-CoA (3-Oxodecanoyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0686 mL	5.3432 mL	10.6864 mL
5 mM	0.2137 mL	1.0686 mL	2.1373 mL
10 mM	0.1069 mL	0.5343 mL	1.0686 mL
50 mM	0.0214 mL	0.1069 mL	0.2137 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481