

DLX-105

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DLX-105 is a human monoclonal antibody (mAb) that targets TNFSF2/TNF α . It reduces epidermal thickness and the expression of K16 and Ki67 in the skin. Additionally, DLX-105 downregulates mRNA levels of IL-17, TNF- α , IL-23p19, IL-12p40, and IFN- γ . This compound is useful for psoriasis research.
Targets(IC50)	IFNAR,Interleukin,TNF

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481