

## Anti-CD88/C5AR1 Antibody (G2\_anti-C5aR)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-CD88/C5AR1 Antibody (G2_anti-C5aR) is a human monoclonal antibody (mAb) targeting C5AR1.
Targets(IC50)	Complement System

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481