

## E3 Ligase Ligand-linker Conjugate 181

## Chemical Properties

CAS No. : 2850354-06-0


Formula: C41H60FN5O12S

Molecular Weight: 866.01

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	E3 ligase Ligand-linker Conjugate 181 is an E3 ubiquitin ligase ligand-linker conjugate (E3 ligase Ligand-linker conjugate). The compound is also known as E3 ligase Ligand-linker ConjugateMS115.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1547 mL	5.7736 mL	11.5472 mL
5 mM	0.2309 mL	1.1547 mL	2.3094 mL
10 mM	0.1155 mL	0.5774 mL	1.1547 mL
50 mM	0.0231 mL	0.1155 mL	0.2309 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481