

GBR-1342

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	GBR-1342 is a human bispecific antibody targeting CD38 and CD3. It inhibits the activity of cancer cell lines Daudi, RAJI, RPMI8226, IM-9, and NCI-H929, and is applicable to multiple myeloma research.
Targets(IC50)	CD38

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481