

SSGJ-613

Chemical Properties

CAS No. : 2775298-78-5

Formula:

Molecular Weight:

Storage: Store at low temperature

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	SSGJ-613 is a fully human anti-IL-1 β monoclonal antibody, which is applicable in the study of gouty arthritis.
Targets(IC50)	Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481