

SCT200

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	SCT200 is a fully humanized IgG1 monoclonal antibody targeting EGFR, with a Kd value of 0.08 nM. It can eliminate tumor cells via complement-dependent cytotoxicity (CDC) and antibody-dependent cellular cytotoxicity (ADCC) through its Fc fragment. SCT200 is utilized in research on refractory RAS and BRAF wild-type metastatic colorectal cancer.
Targets(IC50)	EGFR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481